Trilor FAQ

Like you,
we're experts.

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Bonding questions

Like you,
we're experts.

Empowering
dental professionals
ISO 13845 CERTIFIED | AUTHORISED PARTNER

FAQ

Q Is the bonding step necessary when bonding TriLor to composite?

A Yes, every bit helps.

Q What is the bonding process for TriLor?

A Here is a run-down of the process after milling Trilor:

- Sandblast
- Steam
- Air dry
- Make your restorations using your conventional methods (zir, lithi)
- Now you're ready to put them together!
- Add dip wax to your glazed restoration in the areas that you're not going to have cement
- Apply ceramic primer to the restoration
- Apply composite primer to the TriLor
- Air dry
- Cement crown—use the right cement color!
- Clean up excess cement
- Sandblast
- Add composite primer to TriLor areas
- Add ceramic primer around areas that ceramic will touch the TriLor
- Apply composite
- Cure
- Air Barrier around composite work
- Final cure
- Contour, etc
- Stain glaze
- Light Cure

Q What bonding material can I use on TriLor?

A Some good Resin Cements I and others have used are RelyX by 3M, Multilink by Bisco, Nexus by Kerr and Panavia by Kuraray.

Whatever bonding cement you currently use will work, just be sure it doesn't have any geometric changes like shrinkage or expansion that you hadn't accounted for.

Apply the acrylic resin or composite directly to the TriLor, performing the operative protocol as dictated by the manufacturer of the aesthetic material.

Adhesion protocols for the TriLor is as follows:

1. **Acrylic** | Sandblast the TriLor with Al_2O_3 at 2 Bar (50–110 micron); Air dry well (Oil and water free line); then wet the TriLor with acrylic resin and bond.
2. **Composite** | Sandblast the TriLor with Al_2O_3 at 2 Bar (50–10 micron); Air dry well (Oil and water free line); Prime with composite primer and proceed to apply as you would. Please follow instruction for the composite you use.
3. **Lithium Disilicate** | Sandblast the TriLor with Al_2O_3 at 2 Bar (50–110 micron); Air dry well (Oil and water free line); Sandblast the internal part of the crowns, use etching gel (as recommended) and bond with Universal Resin Cement
4. **Ti-Bases** | Sandblast the TriLor with Al_2O_3 at 2 Bar (50–110 micron); Air dry well (Oil and water free line); Sandblast the intimate part of the Ti-Base, and bond with Universal Resin Cement.

Bond and Cement Details

Esthetic Opportunities

Zirconia, Lithium Disilicate, PMMA, Nano and Composite.

For best results, adhere to material manufacturers supported bonding protocol.

TriLor Treatment

Sandblast TriLor contact areas with Al_2O_3 110 μm /2 bar.

Clean surface with gentle pressurized steam.

All gingival contact areas **should** be sealed (GC OptiGlaze).

Zirconia Bonding

Sandblast intaglio surface with Al_2O_3 110 μm /2 bar.

Clean surface with gentle pressurized steam.

Apply Zirconia Primer (mfr. recommended salinate, i.e. MonoBond Plus, Ivoclar)

Bond with Dual Cure Cement following manufacturer protocols (mfr. recommended, i.e. MultiLink, Ivoclar).

Lithium Disilicate Bonding

Sandblast intaglio surface with Al_2O_3 110 $\mu\text{m}/2$ bar.

Clean surface with gentle pressurized steam.

Etch intaglio surface, Hydrofluoric Acid 5%, 20 seconds.

Rinse for 20 seconds.

Apply Bonding Primer Silane (mfr. recommended, i.e. MonoBond Plus, Ivoclar).

Bond with Dual Cure Cement (mfr. recommended, i.e. Multilink Ivoclar).

PMMA – Milled Bonding

Sandblast intaglio surface with Al_2O_3 110 $\mu\text{m}/2$ bar.

Clean surface with gentle pressurized steam.

Apply primer silane (mfr. recommended, i.e. MonoBond Plus, Ivoclar).

Bond with Dual Cure Cement (mfr. recommended, i.e. Multilink Ivoclar).

Composite/ Nano-Ceramic – Milled Bonding

Sandblast intaglio surface with Al_2O_3 110 $\mu\text{m}/2$ bar.

Clean surface with gentle pressurized steam.

Apply primer silane (mfr. recommended, i.e. MonoBond Plus, Ivoclar).

Bond with Dual Cure Cement (mfr. recommended, i.e. Multilink Ivoclar).

Composite –Build-up Bonding

Sandblast intaglio surface with Al_2O_3 110 $\mu\text{m}/2$ bar.

Clean surface with gentle pressurized steam.

Apply Bonding Adhesive (mfr. recommended, i.e. SR Connect, Ivoclar).

Build-up Composite directly to TriLor.

Q Can TriLor be un-boded from a Ti Base?

A Unfortunately, there is no way of removing the Ti Bases using heat or chemicals. It needs to be physically removed and then rebonded. The trick though is to remove it in a way that it does not compromise the restoration and if it is, it will necessitate a remill of the bridge. The challenge is that if the mistake is during clinical procedures, unfortunately there is no way of fixing it in the laboratory and out of the mouth.

Like you,
we're experts.

Empowering
dental professionals

ISO 13845 CERTIFIED | AUTHORISED PARTNER